




The
CATALYS[®]
Precision Laser System



 has gone **Mobile.**



The time to adopt Laser Cataract Surgery is **NOW.**

Laser Cataract Surgery when you need it.

Deliver excellent patient outcomes through:

- > Gentle docking with minimal IOP rise
- > 3D OCT that guides laser delivery
- > Precise, free floating capsulotomies within 30 μm
- > Complete segmentation and softening of cataract



Partner with ForTec Medical.

On-Sight Mobile Services supports your growth with:

- > No capital outlay or maintenance costs
- > Commitment-free certification
- > Trained and credentialed technicians
- > Professional marketing support
- > Profitable from procedure #1

