



Medical Device ACCESSORIES


The ForTec Smartscope Holmium Laser Fiber is compatible with all open port Holmium lasers. It features a rounded tip which ensures it passes smoothly through the working channel of the ureteroscope.

The fibers are available in 272, 365, 550 and 1000 μm sizes and conveniently packaged sterile for single-use.

- **Smart Design:** The smooth rounded tip passes effortlessly through the scope's working channel reducing damage and costly scope repairs compared to standard flat tip fibers.
- **Improved Performance:** Energy hardens the proprietary rounded tip resulting in a highly durable tip with more consistent energy transmission.
- **Superior/Better Value:** With improved performance, the Smartscope fibers are available at lower price points than standard flat tip fibers providing greater value.



A competitor's traditional flat tip fiber



A SmartScope, Rounded Tip Fiber

[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service