



Case-by-Case
MOBILE SERVICE

UltraPulse CO₂ Laser System

The UltraPulse CO₂ can be used for facial resurfacing for rhytids and acne scarring. The UltraPulse is also used for the removal of perioral and periorbital wrinkles, as well as blepharoplasty.

The UltraPulse CO₂ comes with the CPG (computer Pattern Generator). This helps treat large areas quickly with uniform fluence, allowing precise removal of tissue, layer by layer, with minimal thermal damage.



Laser Specifications:

Laser type: CO₂, Sealed cavity, RF excited, Tankless

Wavelength: 10.6 micron

Power: 1-100W

Dimensions: 77" (h) x 17" (w) x 33" (d)

[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

www.fortecmedical.com