

Imperium


CO₂ Laser Fiber



Medical Device ACCESSORIES

The ForTec Imperium CO₂ Laser Fiber is compatible with all open port CO₂ lasers. It is available in 300 and 500 micron sizes, and packaged sterile for single-use.

- Efficient: 80% power delivery to fiber tip
- Effective: Treatment with high power capacity up to 30W
- Minimal damage: Very limited tissue thermal spread with CO₂ energy
- Ease Of Use: No high pressure gas coolant needed
- Convenient: Renewable fiber tip during surgery



	300μ OTO	300μ	500μ
PSI Rating	30 psi	30 psi	30 psi
Maximum Wattage	10 Watts	30 Watts	30 Watts
Beam Spot Size	200 - 250μ	230μm	320μm
Transmission	60%	80%	80%
Outer Tip Diameter	635μ	740μ	840μ
Aiming Beam	X	X	✓

[THE FORTEC ADVANTAGE]

- ✓ Service and technical expertise
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101