

# Imperium


CO<sub>2</sub> Laser Fiber



## Medical Device ACCESSORIES

The ForTec Imperium CO<sub>2</sub> Laser Fiber is compatible with all open port CO<sub>2</sub> lasers. It is available in 300 and 500 micron sizes, and packaged sterile for single-use.

- Efficient: 80% power delivery to fiber tip
- Effective: Treatment with high power capacity up to 30W
- Minimal damage: Very limited tissue thermal spread with CO<sub>2</sub> energy
- Ease Of Use: No high pressure gas coolant needed
- Convenient: Renewable fiber tip during surgery



## [ THE FORTEC ADVANTAGE ]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

**800.963.7101**

to start ForTec Mobile Service