



Case-by-Case
MOBILE SERVICE

TruBlue

Laser System

The TruBlue laser is a compact and lightweight laser emitting energy up to 10W at 445 nm with peak absorption in hemoglobin. This provides improved coagulation compared to the KTP laser, while also providing good cutting qualities.

It can be used with a variety of fiber sizes in contact or non-contact mode for delicate applications in laryngology, otology and phonosurgery.



[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today


800.963.7101

to start ForTec Mobile Service

TruBlue Laser System



- Excellent photo-angiolysis for improved absorption in hemoglobin
- Permeability in water provides good cutting qualities
- Low depth of penetration for minimal collateral damage
- Highly versatile delivering up to 10W of power with 300, 400 and 600 micron fibers



Truly Flexible



Truly Mobile



Truly Versatile



Truly Noiseless

