



Case-by-Case
MOBILE SERVICE

Single Use **Digital Ureteroscope**

The Pusen flexible ureteroscope delivers high resolution imagery for detailed visualization and accurate navigation for effective treatment of stones.

The single use scope is individually packaged and pre-sterilized to eliminate scope reprocessing, maintenance, repair costs and any compromise in performance which are all common over the life of a reusable scope.

The scope is available in two diameters. The standard scope is 9.5Fr providing a field of view of 120°. A reduced diameter scope of 7.5Fr has a field of view of 90° and provides up to 147% better irrigation



[THE FORTEC ADVANTAGE]

- ✓ Avoid capital and maintenance costs
- ✓ Experienced, professional technician support
- ✓ A stable, professional partner
- ✓ Large inventory of surgical technologies

Call today

800.963.7101

to start ForTec Mobile Service

Single Use Digital Ureteroscope

Specifications:

Deflection of distal tip: 270° up, 270° down


Diameter of insertion position: 7.5 or 9.5Fr

Field of View: 90° (7.5Fr) or 120° (9.5Fr)


Working length: 650mm

Working channel: 3.6Fr


Fully Flexible Design




Lightweight Handle with Deflection Lock




270° Angulation



Slim Design



Bending Section



Irrigation and Laser Locker

